

---

## Product Datasheet

### Human PDGF R alpha / PDGFRA Protein (orb570289)

**Description**

Human PDGF R alpha, His Tag is

(HEK293). It contains AA Gln 24 - Ala ...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

PDGF R alpha

**Preservatives**

PBS, pH7.4

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.4 kDa. The protein migrates as 66-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

NP\_006197.1

**Purity**

95%

**Source**

Human PDGF R alpha, His Tag (orb570289) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Ala 528 (Accession # P16234-1).

**MW**

58.4 kDa

**Uniprot ID**
[P16234](#)
**NCBI**
[NP\\_006197.1](#)
**Expiration Date**

6 months from date of receipt.