

Product Datasheet

Human Tau-441 / 2N4R Protein (orb570300)

Description	Human Tau-441, His Tag is expressed from E. coli cells. It contains AA Met 1 - Leu 441 (Accession # ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Tau
Preservatives	50 mM Tris, 150 mM NaCl, 1 mM TCEP, 1 mM EDTA, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 47.9 kDa. The protein migrates as 60 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_005901.2
Purity	90%
Source	Human Tau-441, His Tag (orb570300) is expressed from E. coli cells. It contains AA Met 1 - Leu 441 (Accession # P10636-8).
MW	47.9 kDa
Uniprot ID	P10636
NCBI	NP_005901.2
Expiration Date	6 months from date of receipt.