


Product Datasheet

Human Thrombospondin-2 Protein (orb570301)

Description	Human Thrombospondin-2, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 19 - ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Thrombospondin-2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 129.7 kDa. The protein migrates as 180-200 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_003238.2
Purity	90%
Source	Human Thrombospondin-2, His Tag (orb570301) is expressed from human 293 cells (HEK293). It contains AA Gly 19 - Ile 1172 (Accession # P35442-1).
MW	129.7 kDa
Uniprot ID	P35442
NCBI	NP_003238.2
Expiration Date	6 months from date of receipt.



SDS-Page analysis of Human Thrombospondi...