

Product Datasheet

Histone H3 (Phospho-Thr3) antibody (orb570424)

Description	Histone H3 (Phospho-Thr3) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HIST1H3A/HIST2H3A/H3F3A
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431, P84243, Q71DI3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.