

Product Datasheet

Histone H3 (Phospho-Thr45) antibody (orb570431)

| | |
|----------------------------|---|
| Description | Histone H3 (Phospho-Thr45) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Synthetic Peptide |
| Target | HIST1H3A/HIST2H3A/H3F3A |
| Preservatives | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P68431 , P84243 , Q71DI3 |
| Dilution Range | WB: 1:1000 |
| Expiration Date | 12 months from date of receipt. |