

Product Datasheet

Histone H1 (Phospho-Thr17) antibody (orb570440)

Description	Histone H1 (Phospho-Thr17) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16403 , P16402 , P10412 , Q02539
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.