

Product Datasheet

TP53 (Phospho-S6) antibody (orb570517)

Description	TP53 (Phospho-S6) antibody
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Synthesized peptide derived from Human p53 around the phosphorylation site of S6.
Target	TP53
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04637
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, ELISA: 1:10000
Expiration Date	12 months from date of receipt.