

## Product Datasheet

### CALM1 (Phospho-T80/S82) antibody (orb570588)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to CALM1 (Phospho-T80/S82)   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB  |
| <b>Immunogen</b>           | Synthesized peptide derived from Human Calmodulin around the phosphorylation site of T80/S82. |
| <b>Target</b>              | CALM1   |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.                  |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                 |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <b>P62158</b>   |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:100-1:300, IF: 1:200-1:1000, ELISA: 1:5000                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |