

Product Datasheet

EFNB1 (Phospho-Y324) antibody (orb570653)

| | |
|----------------------------|---|
| Description | EFNB1 (Phospho-Y324) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | Synthesized peptide derived from Human Ephrin-B1/2/3 around the phosphorylation site of Y324. |
| Target | EFNB1/EFNB2/EFNB3 |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q15768 , P52799 , P98172 |
| Dilution Range | WB: 1:500-1:2000, IHC-P: 1:100-1:300, ELISA: 1:20000 |
| Expiration Date | 12 months from date of receipt. |