

Product Datasheet

PGC (Phospho-Ser15) antibody (orb570848)

Description	PGC (Phospho-Ser15) antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PGC
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20142
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.