

Product Datasheet

CHEK1 (Phospho-S345) antibody (orb570887)

Description	CHEK1 (Phospho-S345) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human Chk1 around the phosphorylation site of S345.
Target	CHEK1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14757
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:10000
Expiration Date	12 months from date of receipt.