

Product Datasheet

EIF2AK2 (Phospho-T451) antibody (orb571211)

Description	EIF2AK2 (Phospho-T451) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human PKR around the phosphorylation site of T451.
Target	EIF2AK2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19525
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:20000
Expiration Date	12 months from date of receipt.