

Product Datasheet

YWHAZ (Phospho-T232) antibody (orb571370)

| | |
|----------------------------|---|
| Description | YWHAZ (Phospho-T232) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | Synthesized peptide derived from Human 14-3-3 $\zeta/6$ around the phosphorylation site of T232. |
| Target | YWHAZ |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P63104 |
| Dilution Range | WB: 1:500-1:2000, IHC-P: 1:100-1:300, ELISA: 1:20000 |
| Expiration Date | 12 months from date of receipt. |