

Product Datasheet

Histone H3 (Phospho-T11) antibody (orb571385)

Description	Histone H3 (Phospho-T11) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from Human Histone H3 around the phosphorylation site of T11.
Target	HIST1H3A/HIST2H3A/H3F3A
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431 , P84243 , Q71DI3
Dilution Range	IF: 1:200-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.