

Product Datasheet

BCL2 (Phospho-S70) antibody (orb571518)

| | |
|----------------------------|---|
| Description | BCL2 (Phospho-S70) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Monkey |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, IP, WB |
| Immunogen | Synthesized peptide derived from Human Bcl-2 around the phosphorylation site of S70. |
| Target | BCL2 |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P10415 |
| Dilution Range | WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, ELISA: 1:10000 |
| Expiration Date | 12 months from date of receipt. |