

Product Datasheet

CDK1 (Phospho-T14) antibody (orb571525)

Description	CDK1 (Phospho-T14) antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Synthesized peptide derived from Human Cdk1/2/3 around the phosphorylation site of T14.
Target	CDK1/CDK2/CDK3
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24941 , Q00526 , P06493
Dilution Range	WB: 1:500-1:2000, IF: 1:200-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.