

Product Datasheet

STAT3 (Phospho-S727) antibody (orb571570)

Description	STAT3 (Phospho-S727) antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Synthesized peptide derived from Human Stat3 around the phosphorylation site of S727.
Target	STAT3
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P40763
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, ELISA: 1:10000
Expiration Date	12 months from date of receipt.