

Product Datasheet

CAMK2A (Phospho-T305) antibody (orb571573)

Description	CAMK2A (Phospho-T305) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthesized peptide derived from Human CaMKII α / β / δ around the phosphorylation site of T305.
Target	CAMK2A/CAMK2B/CAMK2D
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13554 , Q13557 , Q9UQM7
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, ELISA: 1:5000
Expiration Date	12 months from date of receipt.