

Product Datasheet

PPP1R1B (Phospho-Thr75) antibody (orb571674)

Description	PPP1R1B (Phospho-Thr75) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Peptide sequence around phosphorylation site of threonine 75 (A-Y-T(p)-P-P) derived from Human DARPP-32.
Target	PPP1R1B
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UD71
Dilution Range	IF: 1:100-1:200
Expiration Date	12 months from date of receipt.