

Product Datasheet

ARHGAP35 (Phospho-Tyr1105) antibody (orb571743)

Description	ARHGAP35 (Phospho-Tyr1105) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Tyrosine 1105(N-I-Y(p)-S-V) derived from Human GRF-1.
Target	ARHGAP35
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NRY4
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.