

Product Datasheet

RPS6KB1 (Phospho-Thr229) antibody (orb571753)

Description	RPS6KB1 (Phospho-Thr229) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Threonine 229(N-H-Q(p)-G-H) derived from Human p70 S6 Kinase.
Target	RPS6KB1
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23443
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.