

Product Datasheet

CXCR2 (Phospho-Ser347) antibody (orb571780)

Description	CXCR2 (Phospho-Ser347) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 347(R-P-S(p)-F-V) derived from Human IL-8R beta/CDw128 beta
Target	CXCR2
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25025
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.