

Product Datasheet

SMAD2 (Phospho-Ser250) antibody (orb571827)

Description	SMAD2 (Phospho-Ser250) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of serine250(E-L-S(p)-P-T) derived from Human Smad2 .
Target	SMAD2
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15796
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.