

Product Datasheet

CDK2 (Phospho-Thr14) antibody (orb571889)

Description	CDK2 (Phospho-Thr14) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of threonine14(E-G-T(p)-Y-G) derived from Human CDK1/CDC2 .
Target	CDK2
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24941
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.