

Product Datasheet

HSF1 (Phospho-Thr142) antibody (orb571912)

Description	HSF1 (Phospho-Thr142) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Peptide sequence around phosphorylation site of threonine 142(L-L-T(p)-D-V) derived from Human HSF1 .
Target	HSF1
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00613
Dilution Range	IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.