

Product Datasheet

CD40 (Phospho-Thr254) antibody (orb572004)

Description	CD40 (Phospho-Thr254) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of threonine 254 (Q-E-T(p)-L-H) derived from Human TNFRSF5.
Target	CD40
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25942
Dilution Range	WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.