

Product Datasheet

PAK2 (Phospho-Ser141) antibody (orb572008)

Description	PAK2 (Phospho-Ser141) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of serine 141 (Y-L-S(p)-F-T) derived from Human PAK2.
Target	PAK2
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13177
Dilution Range	WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.