

Product Datasheet

MKNK1 (Phospho-Thr255) antibody (orb572018)

Description	MKNK1 (Phospho-Thr255) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of threonine 255 (L-T-T(p)-P-C) derived from Human MNK1.
Target	MKNK1
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BUB5
Dilution Range	WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.