

Product Datasheet

IKBKB (Phospho-Tyr199) antibody (orb572027)

Description	IKBKB (Phospho-Tyr199) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Peptide sequence around phosphorylation site of tyrosine 199(Q-K-Y(p)-T-V) derived from Human IKK- κ .
Target	IKBKB
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14920
Dilution Range	IF: 1:100-1:200
Expiration Date	12 months from date of receipt.