

Product Datasheet

GATA1 (Phospho-Ser310) antibody (orb572046)

| | |
|----------------------------|--|
| Description | GATA1 (Phospho-Ser310) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, WB |
| Immunogen | Peptide sequence around phosphorylation site of serine 310 (K-A-S(p)-G-K) derived from Human GATA1. |
| Target | GATA1 |
| Preservatives | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Form/Appearance | liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P15976 |
| Dilution Range | WB: 1:500-1:1000, IF: 1:100-1:200 |
| Expiration Date | 12 months from date of receipt. |