

## Product Datasheet

# HDAC4 (Phospho-Ser246/259/220) antibody (orb572316)

<b>Description</b>	HDAC4 (Phospho-Ser246/259/220) antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	Peptide sequence around phosphorylation site of serine 246/259/220 (T-A-S(p)-EP) derived from Human HDAC4/HDAC5/HDAC9.
<b>Target</b>	HDAC4/HDAC5/HDAC9
<b>Preservatives</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Form/Appearance</b>	liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q9UQL6</a> , <a href="#">Q9UKV0</a> , <a href="#">P56524</a>
<b>Dilution Range</b>	WB: 1:500-1:1000, IHC-P: 1:50-1:100
<b>Expiration Date</b>	12 months from date of receipt.