

Product Datasheet

Di-Methyl-Histone H3 (Lys14) antibody (orb572419)

Description	Di-Methyl-Histone H3 (Lys14) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide
Target	HIST1H3A/HIST2H3A/H3F3A
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431 , P84243 , Q71DI3
Dilution Range	WB: 1:500-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.