

## Product Datasheet

# RPS6KA1 (Phospho-T359+S363) antibody (orb572475)

<b>Description</b>	Recombinant rabbit monoclonal antibody to RPS6KA1 (Phospho-T359+S363)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IP, WB
<b>Immunogen</b>	A synthesized peptide derived from human Phospho-RPS6KA1 (T359+S363)
<b>Target</b>	RPS6KA1
<b>Preservatives</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q15418</a>
<b>Dilution Range</b>	WB: 1:500-1:5000, IF: 1:20-1:200, IP: 1:200-1:1000
<b>Expiration Date</b>	12 months from date of receipt.