

## Product Datasheet

### NBN (Phospho-S343) antibody (orb572484)

<b>Description</b>	Recombinant rabbit monoclonal antibody to NBN (Phospho-S343)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	A synthesized peptide derived from human Phospho-NBN (S343)
<b>Target</b>	NBN
<b>Preservatives</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">O60934</a>
<b>Dilution Range</b>	WB: 1:500-1:5000
<b>Expiration Date</b>	12 months from date of receipt.