

Product Datasheet

CDK1 (Phospho-T14) antibody (orb572520)

Description	Recombinant rabbit monoclonal antibody to CDK1 (Phospho-T14)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthesized peptide derived from human Phospho-CDK1 (T14)
Target	CDK1/CDK2/CDK3
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P24941 , Q00526 , P06493
Dilution Range	WB: 1:500-1:5000
Expiration Date	12 months from date of receipt.