

Product Datasheet

MAP2K1 (Phospho-T292) antibody (orb572524)

| | |
|----------------------------|---|
| Description | Recombinant rabbit monoclonal antibody to MAP2K1 (Phospho-T292) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | A synthesized peptide derived from human Phospho-MAP2K1 (T292) |
| Target | MAP2K1 |
| Preservatives | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Monoclonal |
| Uniprot ID | Q02750 |
| Dilution Range | WB: 1:500-1:5000 |
| Expiration Date | 12 months from date of receipt. |