

Product Datasheet

CTNNB1 (Phospho-S33/S37) antibody (orb572535)

Description	Recombinant rabbit monoclonal antibody to CTNNB1 (Phospho-S33/S37)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	A synthesized peptide derived from human Phospho-CTNNB1 (S33/S37)
Target	CTNNB1
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P35222
Dilution Range	IF: 1:20-1:200
Expiration Date	12 months from date of receipt.