

Product Datasheet

POLR2A (Phospho-S5) antibody (orb572541)

Description	Recombinant rabbit monoclonal antibody to POLR2A (Phospho-S5)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	A synthesized peptide derived from human Phospho-POLR2A (S5)
Target	POLR2A
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P24928
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:50-1:200, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.