

Product Datasheet

PAK4(Phospho-S474/S560/S602) antibody (orb572557)

Description	Recombinant rabbit monoclonal antibody to PAK4(Phospho-S474/S560/S602)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	A synthesized peptide derived from human Phospho-PAK4/PAK5/PAK6 (S474/S560/S602)
Target	PAK4/PAK5/PAK6
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q9NQU5, Q9P286, O96013
Dilution Range	WB: 1:500-1:5000, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.