

Product Datasheet

LY223982 (orb572866)

Description	LY223982 is a potent and specific inhibitor of leukotriene B4 receptor, with an IC50 of 13.2 nM agai...
Note	For research use only
Formula	$C_{30}H_{30}O_7$
MW	502.56
CAS Number	[117423-74-2]
SMILES	<chem>O=C(O)CCC1=CC(C(C2=CC=CC(C(O)=O)=C2)=O)=CC=C1OCCCC/C=C/C3=CC=C(OC)C=C3</chem>
Expiration Date	12 months from date of receipt.