



## Product Datasheet

### DLinDMA (orb572878)

<b>Description</b>	DLinDMA is a key lipid component of stable nucleic acid lipid particles as a benchmark.
<b>Storage</b>	2 years in -20°C Powder, 2 weeks in 4°C in DMSO, 6 months in -80°C in DMSO
<b>Note</b>	For research use only
<b>Formula</b>	$C_{41}H_{77}NO_2$
<b>MW</b>	616.06
<b>CAS Number</b>	871258-12-7
<b>SMILES</b>	<chem>CCCCC/C=C\C/C=C\C\CCCCCCCCOCC(OCCCCCCCC/C=C\C/C=C\C\CCCC)C N(C)C</chem>
<b>Expiration Date</b>	12 months from date of receipt.



Chemical  
Structure  
(871258-12-7)



CAS: 871258-  
12-7