

## Product Datasheet

### TAS-115 mesylate (orb572929)

<b>Description</b>	TAS-115 mesylate is a potent VEGFR and c-Met/HGFR kinase inhibitor with IC50s of 30 and 32 nM for rV...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{28}H_{27}FN_4O_7S_2$
<b>MW</b>	614.66
<b>CAS Number</b>	[1688673-09-7]
<b>SMILES</b>	<chem>CS(=O)(O)=O.O=C(C1=C(OC)C=C2N=CC=C(OC3=CC=C(NC(NC(CC4=CC=CC=C4)=O)=S)C=C3F)C2=C1)NC</chem>
<b>Expiration Date</b>	12 months from date of receipt.