

Product Datasheet

ML-109 (orb572982)

Description	ML-109 is a potent and full thyroid stimulating hormone receptor (TSHR) agonist, with an EC50 of 40 ...
Note	For research use only
Formula	$C_{31}H_{29}N_3O_5$
MW	523.58
CAS Number	[1186649-91-1]
SMILES	<chem>CC(NC1=CC=C(OCC2=CC(C(N3CC4=CC=CC=C4)NC5=C(C(O)=CC=C5)C3=O)=CC=C2OC)C=C1)=O</chem>
Expiration Date	12 months from date of receipt.