

Product Datasheet

E3 Ligase Ligand 3 (orb573151)

Description	E3 ligase Ligand 3 is a ligand of E3 ligase, used in PROTAC technology.
Note	For research use only
Formula	$C_{15}H_{12}N_2O_7$
MW	332.268
CAS Number	[1061605-21-7]
SMILES	<chem>C(O)(=O)COC1=CC=CC2=C1C(=O)N(C1CCC(=O)NC1=O)C2=O</chem>
Expiration Date	12 months from date of receipt.