

## Product Datasheet

### CP 105696 (orb573184)

<b>Description</b>	CP-105696 is a potent and selective Leukotriene B4 Receptor antagonist, with an IC50 of 8.42 nM.
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>28</sub> H <sub>28</sub> O <sub>4</sub>
<b>MW</b>	428.51952
<b>CAS Number</b>	[158081-99-3]
<b>SMILES</b>	<chem>O=C(C1(C2=CC=C3[C@H](O)[C@@H](CC4=CC=C(C5=CC=CC=C5)C=C4)COC3=C2)CCCC1)O</chem>
<b>Expiration Date</b>	12 months from date of receipt.