

Product Datasheet

SR9238 (orb573186)

Description	SR9238 is a potent and selective LXR inverse agonist.
Note	For research use only
Formula	$C_{31}H_{33}NO_7S_2$
MW	595.73
CAS Number	[1416153-62-2]
SMILES	<chem>CCOC(C1=CC=C(O1)CN(S(=O)(C2=C(C=C(C=C2)C)C)=O)CC3=CC=C(C4=CC(S(C)(=O)=O)=CC=C4)C=C3)=O</chem>
Expiration Date	12 months from date of receipt.