

## Product Datasheet

### GSK2643943A (orb573187)

<b>Description</b>	GSK2643943A is a potent USP20 inhibitor with IC50 = 160 nM.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{17}H_{12}FN_3$
<b>MW</b>	277.3
<b>SMILES</b>	<chem>N#CC1=C(N)NC2=C1C=CC(/C=C/C3=CC=CC(F)=C3)=C2</chem>
<b>Expiration Date</b>	12 months from date of receipt.