

## Product Datasheet

### UT-155 (orb573218)

<b>Description</b>	UT-155 is a selective androgen receptor (AR) antagonist, with a $K_i$ of 267 nM for AR-RBD.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{20}H_{15}F_4N_3O_2$
<b>MW</b>	405.35
<b>CAS Number</b>	[2031161-35-8]
<b>SMILES</b>	<chem>O=C([C@@](C)(O)CN1C(C=CC(F)=C2)=C2C=C1)NC3=CC=C(C#N)C(C(F)(F)F)=C3</chem>
<b>Expiration Date</b>	12 months from date of receipt.