

Product Datasheet

IWP-O1 (orb573386)

Description	IWP-O1 is a porcn inhibitor with an EC50 value of 80 pM in a cultured cell reporter assay of Wnt sig...
Note	For research use only
Formula	$C_{26}H_{20}N_6O$
MW	432.48
CAS Number	[2074607-48-8]
SMILES	<chem>C1=CC=C(C=C1)C2=CN=C(C=C2)NC(=O)CN3C(=C(N=N3)C4=CC=NC=C4)C5=CC=CC=C5</chem>
Expiration Date	12 months from date of receipt.