

## Product Datasheet

### NSC617145 (WRN inhibitor) (orb573410)

|                        |   |
|------------------------|---|
| <b>Description</b>     | NSC617145 is a Werner syndrome helicase (WRN) helicase inhibitor (IC50 = 250 nM). |
| <b>Note</b>            | For research use only   |
| <b>Formula</b>         | $C_{13}H_{10}Cl_4N_2O_4$  |
| <b>MW</b>              | 400.04  |
| <b>CAS Number</b>      | [203115-63-3]   |
| <b>SMILES</b>          | <chem>CC(CN1C(C(Cl)=C(Cl)C1=O)=O)(C)CN2C(C(Cl)=C(Cl)C2=O)=O</chem>                |
| <b>Expiration Date</b> | 12 months from date of receipt.   |